

# USP16 Antibody



PACO13121

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

USP16

**Uniprot**

Q9Y5T5

**Synonyms:**

ubiquitin specific peptidase 16; USP16; UBP-M ;

**Immunogen:**

Human USP16.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A