

USP2 Antibody



PACO13122

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

USP2

Uniprot

O75604

Synonyms:

ubiquitin specific peptidase 2; USP2; UBP41; USP9 ;

Immunogen:

Human USP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A