

USP22 Antibody



PACO13124

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

USP22

Uniprot

Q9UPT9

Synonyms:

ubiquitin specific peptidase 22; USP22; KIAA1063; USP3L ;

Immunogen:

Human USP22.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A