

# USP36 Antibody



PACO13132

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

USP36

**Uniprot**

Q9P275

**Synonyms:**

ubiquitin specific peptidase 36; USP36; DUB1 ;

**Immunogen:**

Human USP36.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A