

USP45 Antibody



PACO13137

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

USP45

Uniprot

Q70EL2

Synonyms:

ubiquitin specific peptidase 45; USP45; MGC14793 ;

Immunogen:

Human USP45.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A