

# USP46 Antibody



PACO13138

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

USP46

**Uniprot**

P62068

**Synonyms:**

ubiquitin specific peptidase 46; USP46; FLJ11850; FLJ12552; FLJ14283; FLJ39393 ;

**Immunogen:**

Human USP46.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A