

USP49 Antibody



PACO13140

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

USP49

Uniprot

Q70CQ1

Synonyms:

ubiquitin specific peptidase 49; USP49; MGC20741 ; Deubiquitinating enzyme 49; Ubiquitin thiolesterase 49

Immunogen:

Human USP49.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A