

UTP23 Antibody



PACO13148

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UTP23

Uniprot

Q9BRU9

Synonyms:

UTP23, small subunit (SSU) processome component, homolog (yeast); UTP23; C8orf53; MGC14595 ;

Immunogen:

Human UTP23.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A