

VAC14 Antibody



PACO13153

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VAC14

Uniprot

Q08AM6

Synonyms:

Vac14 homolog (S. cerevisiae); VAC14; ArPIKfyve; FLJ10305; FLJ36622; FLJ46582; MGC149815; MGC149816; TAX1BP2; TRX ;

Immunogen:

Human VAC14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A