

# VANGL1 Antibody



PACO13159

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VANGL1

**Uniprot**

Q8TAA9

**Synonyms:**

vang-like 1 (van gogh, Drosophila); VANGL1; LPP2; MGC5338; STB2; STBM2 ;

**Immunogen:**

Human VANGL1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A