

PACO13162

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VARs

Uniprot

P26640

Synonyms:

valyl-tRNA synthetase; VARs; G7A; VARs1; VARs2 ;

Immunogen:

Human VARs.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A