

# VASH2 Antibody



PACO13164

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VASH2

**Uniprot**

Q86V25

**Synonyms:**

vasohibin 2; VASH2; FLJ12505; RP11-275G3.1 ;

**Immunogen:**

Human VASH2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A