

VAX1 Antibody



PACO13167

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

VAX1

Uniprot

Q55QQ9

Synonyms:

ventral anterior homeobox 1; VAX1; MGC126743; MGC126745 ;

Immunogen:

Human VAX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A