

VCP Antibody



PACO13169

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VCP

Uniprot

P55072

Synonyms:

valosin-containing protein; VCP; IBMPFD; MGC131997; MGC148092; MGC8560; TERA; p97 ;

Immunogen:

Human VCP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A