

# VEGFA Antibody



PACO13174

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VEGFA

**Uniprot**

P15692

**Synonyms:**

vascular endothelial growth factor A; VEGFA; MGC70609; MVCD1; VEGF; VEGF-A; VPF ;

**Immunogen:**

Human VEGFA.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A