

VEZT Antibody



PACO13175

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VEZT

Uniprot

Q9HBM0

Synonyms:

vezatin, adherens junctions transmembrane protein; VEZT; DKFZp761C241; VEZATIN ;

Immunogen:

Human VEZT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A