

VMO1 Antibody



PACO13180

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VMO1

Uniprot

Q7Z5L0

Synonyms:

vitelline membrane outer layer 1 homolog (chicken); VMO1; ERGA6350; MGC125880; MGC125881; PRO21055 ;

Immunogen:

Human VMO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A