

VPS41 Antibody



PACO13195

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VPS41

Uniprot

P49754

Synonyms:

vacuolar protein sorting 41 homolog (S. cerevisiae); VPS41; HVPS41; HVSP41; hVps41p ;

Immunogen:

Human VPS41.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A