

VPS72 Antibody



PACO13201

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VPS72

Uniprot

Q15906

Synonyms:

vacuolar protein sorting 72 homolog (S. cerevisiae); VPS72; CFL1; Swc2; TCFL1; YL-1; YL1 ;

Immunogen:

Human VPS72.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A