

VSNL1 Antibody



PACO13206

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VSNL1

Uniprot

P62760

Synonyms:

visinin-like 1; VSNL1; HLP3; HPCAL3; HUVISL1; VILIP; VILIP-1 ;

Immunogen:

Human VSNL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A