

VWA1 Antibody



PACO13211

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VWA1

Uniprot

Q6PCB0

Synonyms:

von Willebrand factor A domain containing 1; VWA1; DKFZp761O051; FLJ22215; VWA-1; WARP ;

Immunogen:

Human VWA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A