

# VWA2 Antibody



PACO13212

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VWA2

**Uniprot**

Q5GFL6

**Synonyms:**

von Willebrand factor A domain containing 2; VWA2; AMACO; CCSP-2; MGC157755 ;

**Immunogen:**

Human VWA2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A