

PACO13218

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

WASL

Uniprot

O00401

Synonyms:

Wiskott-Aldrich syndrome-like; WASL; DKFZp779G0847; MGC48327; N-WASP; NWASP ;

Immunogen:

Human WASL.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A