

WDFY3 Antibody



PACO13223

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WDFY3

Uniprot

Q8IZQ1

Synonyms:

WD repeat and FYVE domain containing 3; WDFY3; ALFY; KIAA0993; MGC16461; ZFYVE25 ;

Immunogen:

Human WDFY3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A