

PACO13224

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WDFY4

Uniprot

Q6ZS81

Synonyms:

WDFY family member 4; WDFY4; C10orf64; FLJ17416; FLJ36288; FLJ45748; KIAA1607; MGC40604; MGC45899 ;

Immunogen:

Human WDFY4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A