

# WDR18 Antibody



PACO13227

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

WDR18

**Uniprot**

Q9BV38

**Synonyms:**

WD repeat domain 18; WDR18; MGC2436; R32184\_1 ;

**Immunogen:**

Human WDR18.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A