

WDR24 Antibody



PACO13230

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

WDR24

Uniprot

Q96S15

Synonyms:

WD repeat domain 24; WDR24; C16orf21; DKFZp434F054; JFP7 ;

Immunogen:

Human WDR24.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A