

WDR5 Antibody



PACO13236

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

WDR5

Uniprot

P61964

Synonyms:

WD repeat domain 5; WDR5; BIG-3; SWD3 ;

Immunogen:

Human WDR5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A