

WDR6 Antibody



PACO13237

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

WDR6

Uniprot

Q9NNW5

Synonyms:

WD repeat domain 6; WDR6; FLJ10218; FLJ52552; FLJ56107; MGC126756; MGC142027 ;

Immunogen:

Human WDR6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A