

WHSC1L1 Antibody



PACO13241

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NSD3

Uniprot

Q9BZ95

Synonyms:

Wolf-Hirschhorn syndrome candidate 1-like 1; WHSC1L1; DKFZp667H044; FLJ20353; MGC126766; MGC142029; NSD3; pp14328 ;

Immunogen:

Human WHSC1L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A