

PACO13243

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

WIPI2

**Uniprot**

Q9Y4P8

**Synonyms:**

WD repeat domain, phosphoinositide interacting 2; WIPI2; Atg21; CGI-50; DKFZp434J154; DKFZp686P02188; FLJ12979; FLJ14217; FLJ42984; WIPI-2 ;

**Immunogen:**

Human WIPI2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A