

WISP1 Antibody



PACO13244

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WISP1

Uniprot

O95388

Synonyms:

WNT1 inducible signaling pathway protein 1; WISP1; CCN4; WISP1c; WISP1i; WISP1tc ;

Immunogen:

Human WISP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A