

WNT2 Antibody



PACO13245

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

WNT2

Uniprot

P09544

Synonyms:

wingless-type MMTV integration site family member 2; WNT2; INT1L1; IRP ;

Immunogen:

Human WNT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A