

WSB1 Antibody



PACO13251

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

WSB1

Uniprot

Q9Y6I7

Synonyms:

WD repeat and SOCS box-containing 1; WSB1; SWIP1; WSB-1 ;

Immunogen:

Human WSB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A