

PACO13255

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

WWOX

**Uniprot**

Q9NZC7

**Synonyms:**

WW domain containing oxidoreductase; WWOX; D16S432E; FOR; FRA16D; HHCMA56; PRO0128; SDR41C1; WOX1 ;

**Immunogen:**

Human WWOX.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A