

WWP2 Antibody



PACO13257

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

WWP2

Uniprot

O00308

Synonyms:

WW domain containing E3 ubiquitin protein ligase 2; WWP2; AIP2; WWp2-like ;

Immunogen:

Human WWP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A