

XPO1 Antibody



PACO13264

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

XPO1

Uniprot

O14980

Synonyms:

exportin 1 (CRM1 homolog, yeast); XPO1; CRM1; DKFZp686B1823 ;

Immunogen:

Human XPO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A