

XPO6 Antibody



PACO13265

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

XPO6

Uniprot

Q96QU8

Synonyms:

exportin 6; XPO6; EXP6; FLJ22519; KIAA0370; RANBP20 ;

Immunogen:

Human XPO6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A