

XRCC4 Antibody



PACO13268

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

XRCC4

Uniprot

Q13426

Synonyms:

X-ray repair complementing defective repair in Chinese hamster cells 4;

Immunogen:

Human XRCC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A