

YARS2 Antibody



PACO13272

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

YARS2

Uniprot

Q9Y2Z4

Synonyms:

tyrosyl-tRNA synthetase 2, mitochondrial; YARS2; CGI-04; FLJ13995; TYRRS; mt-TyrRS ;

Immunogen:

Human YARS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A