

# YEATS4 Antibody



PACO13275

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

YEATS4

**Uniprot**

O95619

**Synonyms:**

YEATS domain containing 4; YEATS4;4930573H17Rik; B230215M10Rik; GAS41; NUBI-1; YAF9 ;

**Immunogen:**

Human YEATS4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A