

YPEL5 Antibody



PACO13281

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

YPEL5

Uniprot

P62699

Synonyms:

yippee-like 5 (Drosophila); YPEL5; CGI-127 ;

Immunogen:

Human YPEL5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A