

YTHDF1 Antibody



PACO13282

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

YTHDF1

Uniprot

Q9BYJ9

Synonyms:

YTH domain family, member 1; YTHDF1; C20orf21; FLJ20391 ;

Immunogen:

Human YTHDF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A