

# ZBTB40 Antibody



PACO13291

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZBTB40

**Uniprot**

Q9NUA8

**Synonyms:**

zinc finger and BTB domain containing 40; ZBTB40; KIAA0478; MGC133098 ;

**Immunogen:**

Human ZBTB40.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A