

ZC3HC1 Antibody



PACO13295

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZC3HC1

Uniprot

Q86WB0

Synonyms:

zinc finger, C3HC-type containing 1; ZC3HC1;1110054L24Rik; HSPC216; NIPA ;

Immunogen:

Human ZC3HC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A