

ZDHHC17 Antibody



PACO13297

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ZDHHC17

Uniprot

Q8IUH5

Synonyms:

zinc finger, DHHC-type containing 17; ZDHHC17; HIP14; HIP3; HSPC294; HYPH; KIAA0946 ;

Immunogen:

Human ZDHHC17.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A