

ZFP36L2 Antibody



PACO13301

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZFP36L2

Uniprot

P47974

Synonyms:

zinc finger protein 36, C3H type-like 2; ZFP36L2; BRF2; ERF-2; ERF2; RNF162C; TIS11D ;

Immunogen:

Human ZFP36L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A