

ZG16 Antibody



PACO13306

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ZG16

Uniprot

O60844

Synonyms:

zymogen granule protein 16 homolog (rat); ZG16; FLJ43571; FLJ92276; JCLN; JCLN1; MGC183567; MGC34820; ZG16A; hZG16 ;

Immunogen:

Human ZG16.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A