

# ZMAT3 Antibody



PACO13310

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZMAT3

**Uniprot**

Q9HA38

**Synonyms:**

zinc finger, matrin type 3; ZMAT3; FLJ12296; MGC10613; PAG608; WIG-1; WIG1 ;

**Immunogen:**

Human ZMAT3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A