

ZKSCAN8 Antibody



PACO13318

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZKSCAN8

Uniprot

Q15776

Synonyms:

zinc finger protein 192; ZNF192; LD5-1; ZKSCAN8 ;

Immunogen:

Human ZNF192.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A